Auto SureFloat™

**Floating Pressure Management
with Alternating Pressure and Active Sensor Technology™**

- REF C1000MS Control Unit
- REF M1000S Series Mattress

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- REF C1000MS regeleenheid
- Matras REF M1000S Series

Table Of Contents

<u>Section</u>	<u>Description</u>	<u>Page</u>
1.0	Warranty.....	2
2.0	Symbols.....	2
3.0	Indications for Use.....	2
4.0	Contraindications.....	2
5.0	Safety Precautions.....	3
6.0	Description.....	4
7.0	Features	4
8.0	Operator Control Panel.....	6
9.0	Setting Up the System.....	7
10.0	Instructions for Use	7
11.0	Transferring Patients In and out of Bed.....	9
12.0	CPR (Cardio Pulmonary Resuscitation).....	9
13.0	System Cleaning, Care and Storage.....	10
14.0	Troubleshooting.....	11
15.0	Specifications, Control Unit.....	12

Before you begin... **Important**

Before using the Auto SureFloat™ System, please read and understand this manual and all safety precautions prior to each application.

Only qualified medical service personnel should attempt to repair this device. A *Service Manual* is available for use by qualified personnel. For assistance contact your local dealer. If additional assistance is needed, contact Gaymar's technical service department.

International Customers:

Direct (716) 662-8636
Fax: (716) 662-0730

Domestic Customers:

Toll Free (800) 828-7341 extension 739
Fax: (716) 662-8795

1.0 Warranty

The Control Unit is warranted free of defects in material and workmanship for a period of two (2) years.

The Mattress is warranted free of defects in material and workmanship each for a period of one (1) year.

The Control Unit, and Mattress are warranted under the terms and conditions of the Gaymar warranty in place at the time of purchase. A copy of the warranty is available upon request. Gaymar disclaims all implied warranties including, but not limited to, the implied warranties of merchantability and of fitness for a particular purpose.

Control Units may be returned to the factory for servicing. Contact Gaymar Technical Service for return authorization prior to return.

Please contact Gaymar Customer Service if you have warranty questions.

Direct (716) 662-2551 extension 739
Toll Free (800) 828-7341

2.0 Symbols

Attention, consult accompanying documents

Type BF equipment

Dangerous voltage

Protective earth

3.0 Indications for Use

This device is intended to assist in treating and preventing pressure ulcers.

4.0 Contraindications

Air support treatment is not recommended when spinal stability is a concern.

5.0 Safety Precautions

Review the following SAFETY PRECAUTIONS prior to using the Auto SureFloat™ System.

DANGER

Risk of electric shock. Refer servicing to qualified service personnel.

WARNING

- Disinfect the Auto SureFloat™ system between patient installations. Failure to disinfect may risk cross-contamination and infection.
- Check patient at least every 8 hours or once per shift, to assure proper system operation.
- Deflate mattress before performing CPR or CPR will be ineffective.

CAUTION

- Some medical conditions may not respond to treatment of this type. Patient's skin condition should be inspected regularly. Consult physician if any redness or skin breakdown occurs.
- For grounding reliability, plug only into a properly grounded outlet.
- Make certain all mattress straps are secured to the bed frame to prevent mattress from sliding and causing patient injury.
- Raise the bed side rails before initiating treatment to prevent possible injury to patient.
- Use minimal layers of sheeting and incontinence pads. Too many layers between the patient's skin and the support surface will reduce the effectiveness of the system.
- Do not pull linens tightly over mattress. Tight sheets can cause "hammocking" and reduce effectiveness of treatment.
- Medical Electrical Equipment requires special precautions regarding EMC and should be installed and put into service according to the EMC information provided in the accompanying documents.
- Portable and RF communications equipment can affect medical electrical equipment.

6.0 Description

The Auto SureFloat™ System is a portable air flotation system with on-demand alternating low pressure treatment, designed to provide benefit to patients suffering from, or at risk of, developing pressure ulcers. The system, consisting of a Control Unit and Air Flotation Mattress with it's Top Sheet, is designed to provide pressure management and patient comfort.

The Control Unit inflates and maintains the Mattress to user-selected settings, and can also provide on demand alternating low pressure treatment if desired. Used with the M1000S mattress, featuring Active Sensor Technology, the system automatically adjusts the Mattress to the optimal fill level ensuring complete flotation pressure management.

Mattress

The Mattress base consists of twenty transverse air cells constructed of a low shear nylon. Each of the individual cells is a minimum of eight inches tall (inflated height) and spatially oriented above a two inch convoluted foam base. The M1000S Mattress utilizes 4 specialized air cells in the center region of the Mattress which contain active sensors for detecting patient flotation level.

Covering the entire Mattress assembly is a low-friction, low-shear producing, vapor permeable, nylon Top Sheet. The Top Sheet surface creates a membrane impermeable to liquids, air and bacteria though still permeable to water vapor. Quilted onto the base of the Mattress Top Sheet is a spun bonded polyester fiberfill that exhibits less frictional resistance to nylon than that of the patient's skin. When the patient moves, the Top Sheet tends to move relative to the air cells rather than relative to the patient, minimizing shearing effects. Additionally, the quilting acts to provide a diffusion layer by which the water vapor can be drawn away from the patient, resulting in more efficient evaporation and subsequent reduced skin maceration.

7.0 Features


7.1 C1000MS Control Unit

- Max inflate for easy patient transfer.
- Soft/firm comfort control adjustment for customized treatment.
- On-demand, alternating, low pressure management.
- Automatic patient flotation control via Active Sensor Technology™ (when used with the M1000S mattress) to ensure optimal pressure management.
- CPR quick release providing for rapid deflation.
- Sleek and durable plastic housing is extremely compact and lightweight.
- Non-slip mounting brackets for secure placement on the foot board of bed.
- Illuminated LED control panel for easy viewing.
- Control panel Lockout Mode to prevent inadvertent patient setting changes.
- Alarm Silence feature to temporarily suspend audible alarm.

7.2 M1000 Series Air Flotation Mattress

- Horizontal, tubular cell construction for optimal comfort and pressure management.
- Individual air cells designed for maximum support and pressure distribution.
- High density foam base for additional protection.
- Low friction, low shear Top Sheet to aid in prevention of skin breakdown.
- Dimensions - base model 10" x 34" x 80" (25 cm x 86 cm x 203 cm) when fully inflated. Other sizes available.
- M1000S model air cells contain active sensors for automatic patient flotation control.
- Optional Footstop/ZIP™ cell to minimize heel interface pressure while maintaining proximal positioning of feet.
- Removable cells for greater flexibility and easy maintenance.

8.0 Operator Control Panel

- 1.** **ON/STANDBY** switches Control Unit between On and STANDBY modes.
- 2.** **POWER INDICATOR LED** indicates AC power is present to the Control Unit.
- 3.** **MODE SELECT** Press repeatedly to toggle through operating modes, Auto/Float, Auto/A/P, Manual/Float and Manual/A/P.

3.1 AUTO MODE (Sensor cord from M1000S mattress must be connected to control unit) Mattress will automatically inflate to optimal level based on information received from active sensors in mattress. The indicator adjacent to the icon, will light while automatic mode is active.

Note: Connecting the mattress sensor cord while the unit is not in automatic mode will cause the Control Unit to enter automatic mode immediately.

If Sensor cord from mattress becomes disconnected from control unit while in Auto mode, alarm will sound.

3.2 MANUAL MODE Mattress inflation level user-determined by Comfort Adjust (- , +) setting. The indicator adjacent to the icon will light while manual mode is active.

3.3 FLOAT MODE all cells within mattress will be uniformly inflated when this icon is illuminated.

3.4 AP MODE. (Alternating Pressure) Adjacent cells within mattress will be alternately inflated and deflated in this mode. The indicator adjacent to the

icon will light while AP mode is active. Note: AP mode is temporarily suspended in max inflate mode.

- 4.** **SOFT (-)** decreases the inflation level of the bed while Control Unit is in manual mode only. Pressing once will decrease the setting by one comfort level. Pressing and holding will rapidly change the setting.
- 5.** **COMFORT CONTROL INDICATOR** displays the current inflation level setpoint.
- 6.** **FIRM (+)** increases inflation level of the bed while in manual mode only. Pressing once will increase the setting by one comfort level setting (~2.5mmHg). Pressing and holding will rapidly change the setting.
- 7.** **MAX INFLATE** selects/deselects the max inflate mode (maximum mattress firmness). The indicator adjacent to the button will flash while max inflate mode is active. After 30 minutes, max inflate mode will cancel.
- 8.** **LOCKOUT** activates or cancels lockout mode. Press and hold for three seconds to activate or cancel mode. The indicator adjacent to the button will light while lockout is active. In lockout mode, only the MAX INFLATE button is operational. All other buttons are temporarily inactive.
- 9.** **ALARM SILENCE** silences indefinitely, the audible alarm.
- 10.** **POWER FAIL INDICATOR** flashes when power supplied from line cord is interrupted while the Control Unit is operating. To avoid an inadvertent power fail alarm, place the unit in STANDBY mode before disconnecting power from Control Unit.

9.0 Setting Up the System

Mattress

1. Place the mattress on the bed frame with the connecting hoses at the foot of the bed.
2. Before securing the mattress straps, raise the head and foot ends of the bed frame. Raise the side rails. Loosely, secure mattress to bed frame by looping the nylon straps on mattress through the bed frame and then back through the D-rings on mattress. Make sure the side rails can still move up and down without interference.
3. Lower frame to its horizontal position and tighten straps.
4. Cover the mattress with the low-shear Top Sheet and attach to the mattress by looping the Velcro straps on the Top Sheet through the D-rings on the mattress. Verify Top Sheet is loosely fit so as not to "hammock" patient above mattress.

Control Unit / Linens / Patient

1. Suspend Control Unit on footboard of bed using the mounting brackets on back of unit.
2. Connect the dual quick-connect couplings of the mattress hoses to the mating connectors located on the side of the Control Unit.
3. Check that the hoses are secure by gently pulling. Ensure the hoses are not kinked or tucked under the mattress.
4. If using an M1000S Mattress, connect the plug end of the mattress sensor cord into the jack on the side of the Control Unit.
5. Plug the AC power cord into properly grounded outlet. The Control Unit will enter STANDBY mode.
6. Press the ON/STANDBY button on the lower right of the control panel. The mattress will begin to inflate. This could take approximately 15 minutes for complete inflation.
7. Apply hospital linens and/or incontinence pad over the top of the Mattress. Linens should be loose to prevent "hammocking".
8. Place patient on Mattress and position along center of bed. Raise the bed side rails.

10.0 Instructions for Use

10.1 Use With M1000 Mattress

Manual Mode – Float

1. Press the Mode Select button until the icons adjacent to the MANUAL and FLOAT icons are both illuminated
2. Adjust the inflation level for patient comfort using the SOFT (-) and FIRM (+) buttons. These buttons adjust the inflation pressures in the mattress from a minimum of 7.5 mmHg (one LED lighted on the COMFORT CONTROL INDICATOR) to a maximum of 30 mmHg (all LED's lighted).
3. Manually ensure the patient is sufficiently supported. In the mattress, place your hand between two air cells and then under the patient's pelvic region. For proper inflation, there should be clearance for 3-4 fingers (stacked vertically) beneath the patient's pelvic region.

WARNING

Check patient at least every 8 hours or once per shift, to assure proper system inflation.

4. To prevent accidental changes in settings, press and hold the LOCKOUT button for approximately 3 seconds.

Manual Mode – A/P

5. If alternating low pressure management is desired, press the Mode Select button until the lights adjacent to the MANUAL and A/P icons are illuminated. If A/P (alternating pressure) is selected in combination with Manual mode, every other cell within the mattress will alternately increase and decrease pressure every 2½ minutes.

10.2 Use with M1000S Mattress (with Active Sensor Technology™)**Automatic Mode - Float**

1. Automatic mode of operation is engaged automatically at the press of the ON button (when the M1000S Mattress sensor cord is plugged into the Control Unit). The indicator adjacent to the AUTO icon will light indicating the automatic mode is active.
2. To select Float, press the Mode Select button until the indicators adjacent to the Float and Auto icons are illuminated.
3. The automatic mode requires **no** user input to determine the appropriate inflation level for the patient. The active sensors in the mattress located under the patient's pelvic region, in conjunction with the Control Unit, automatically determine the optimal inflation level for the patient to ensure complete flotation.

Theory

While in Automatic Float mode, the Control Unit continually repeats a deflation/inflation cycle. The inflation level of the mattress is lowered until the lowest allowable flotation height is sensed (via feedback from the active sensors in the mattress). It is then slowly increased until the appropriate patient flotation height is sensed. This is approximately four fingers width under the patient's pelvic region but will vary with patient weight and build.

If the patient sinks into the mattress for any reason, (such as elevation of the head section or other part of the bed) the Control Unit will automatically increase the inflation level until the appropriate flotation height is again realized.

Disconnecting the mattress sensor cord from the Control Unit while in automatic mode will cause an audible alarm to sound and indicators adjacent to the ALARM SILENCE and AUTO icons to flash. This alarm can be indefinitely silenced by pressing the ALARM SILENCE button. Reconnect the mattress sensor cord to the jack on the right side of the unit, or press the Mode Select button to disengage the automatic mode.

Connecting the mattress sensor cord to the Control Unit while the system is not in automatic mode will engage the automatic mode automatically and light the indicator adjacent to the AUTO button

4. To prevent accidental changes in settings, press and hold the LOCKOUT button for approximately 3 seconds.

Automatic Mode - AP

1. If alternating low pressure therapy is also desired, press the Mode Select button until the indicators adjacent to the Auto and A/P icons are illuminated.
2. In Auto – A/P mode, alternate cells within the mattress are inflated and deflated every 2½ minutes. At the end of each 2½ minute cycle the system will evaluate the floatation level, (determined by the automatic sensors within the mattress) and make comfort level adjustments as necessary before resuming the A/P operation.

11.0 Transferring Patients In and Out of Bed

1. For patient transfer press the MAX INFLATE button. This will maximally inflate the mattress to provide a firm and stable support surface (and suspend AP and automatic modes). Pressing the MAX INFLATE button again will cancel this mode and return the mattress to previously chosen settings.
2. While in the max inflate mode, the indicator adjacent to the MAX INFLATE button will flash to remind the caregiver that the system is in this state. After approximately 30 minutes, the Control Unit will automatically cancel the max inflate mode and return to the previously chosen settings.

12.0 CPR (Cardio Pulmonary Resuscitation)

To deflate mattress for CPR:

1. Disconnect the dual quick-connect couplings of the mattress hose from the Control Unit. Mattress will deflate in approximately 15-20 seconds.
2. Proceed with CPR procedures.

13.0 System Cleaning, Care and Storage

13.1 Cleaning

WARNING

Disconnect the AC power cord from the wall outlet before attempting to clean the Control Unit. Do not heat or steam autoclave any component of the system.

1. To clean, use soap, water and a clean cloth to wipe down the Control Unit, power cord, hoses and mattress. Do not use abrasive cleaners on the mattress. Wipe dry with a clean, dry cloth. Note: Blood and other body fluids must be thoroughly cleaned from all surfaces before applying disinfectants.
2. Apply any EPA approved disinfectant to the external surfaces of the Control Unit, hoses and mattress. Allow contact time according to manufacturers instructions. The solution contact time is what makes disinfection effective.
3. Wipe down the mattress with a clean, dry cloth to remove any excess disinfectant.
4. Top sheets of mattress may also be laundered between patient uses or as required to maintain good patient hygiene. Fill the washing machine with warm water (70 -140 °F or 21 - 60 °C). Add one cup of laundry detergent. Place no more than four top sheets in a single extra large load capacity washing machine. When wash cycle is complete, remove promptly from machine and ensure all excess water is drained from load. Place dryer on LOWEST heat setting, or AIR FLUFF if available until dry. Verify top sheet is completely dry before placing under patient.
5. If individual air cells of mattress become soiled, clean and disinfect as described above or simply replace air cell with a clean replacement. Single air cell replacement can be successfully achieved with patient remaining on the mattress. Air cells which do not contain active sensors can also be laundered. A plug is available from Gaymar to prevent water from getting inside these during laundering (P/N 30287).

WARNING

Do not launder the Claret colored Active Sensor Technology™ air cells of the M1000S mattress. Damage may result. If one or more of the Active Sensor Technology™ air cells become damaged, they can only be replaced as a complete set.

13.2 Storage

Control Unit

1. Check the AC power cord and plug for abrasions or excessive wear.
2. Coil power cord up and place both cord and Control Unit in plastic bag for storage.

Mattress

1. Check the air manifold of the mattress for kinks or breaks. Replace if necessary.
2. To deflate, disconnect hose assembly from Control Unit.
3. Coil the hose assembly loosely and roll the deflated air mattress starting at the head end. Use the straps for containment. Place in plastic bag for storage.

14.0 Troubleshooting

Symptom:	Corrective Action
Mattress does not inflate or is too soft.	Make certain the air hoses are not kinked, cracked, or split. Verify the hose couplings are fully engaged to controller. Verify that all 20 air cells are connected to the air manifold within mattress.
Power loss indicator flashing and audible alarm sounding	Power has been disconnected from unit while in an operational mode. Restore power, put unit in STANDBY mode, then remove power. Unit must always be returned to STANDBY mode before disconnecting power.
Indicator adjacent to ALARM SILENCE button is flashing and audible alarm is sounding.	If in AUTO mode, check the connection of the M1000S mattress sensor cord to the jack on the side of the unit. Reconnect if necessary. Verify hose couplings are fully engaged to controller. Verify that all 20 air cells are connected to the air manifold within mattress.

15.0 Specifications, Control Unit

	C1000MS
Enclosure Dimensions	10" x 12" x 4-1/2" (25 cm x 30 cm x 11 cm)
Weight	9.5 pounds (4.3 kg)
Power Cord	Detachable 14' , #18 AWG minimum, with ground wire
Overcurrent Protection	Primary 2.5A, 250V, T
Input	120 -230 VAC, 50/60Hz, 120V – 0.44A, 230V – 0.33A
Operating Ambient Temperature Range	60 to 85°F (16 to 29°C)
Classification	<p>Class I grounded equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.</p> <p>Type BF equipment</p>  <p>UL2601-1, CSA C22.2 NO. 601.1 EN 60601-1</p> <p>IPXO, enclosed equipment without protection against ingress of water.</p> <p>Continuous operation</p>
Electromagnetic Compatibility	Meets EN60601-1-2:2001 (CISPR 11 Classified as Class B, Group 1 ISM equipment)

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INVATECH GmbH & Co.
Sunmed Medizintechnische Produkte GmbH